

DermaLite


Ryan Caudill, Theo Villalva, & DJ Ravenell

DermaLite

Create an mobile app where a user can take a picture of a mole and be told:


- If it is a mole
- Benign or Malignant
- Type of malignancy

All done on-device for privacy and efficiency




Method

- Original plan
 - Single model, show output to user
 - Inaccurate, unreliable
- New plan




Mole Detector

- Filter out pictures of non-moles
- Challenges
 - Initial Attempt
 - Would not work with any non-dataset image
 - Determined on dermoscopic imaging side-effects
 - Lens distortion, high-quality imaging, scales
 - Second Attempt
 - Heavy preprocessing
 - Cropping, Gaussian blur and noise, grayscale




Mole Detector

- Results seemed quite effective
- Limitations
 - Skin color
 - Size of mole
 - Type of lesion
 - Warts, birthmarks, burns



Benign vs Malignant

- ResNet-50
- Challenges
 - Naturally more difficult and specific
 - Limited to HAM10000 and ISIC 2024
 - 3x benign vs malignant case
 - Malignant is more important
 - Training would plateau
- Solution
 - Preprocessing (same as previous model)
 - Replaced FC layer with random forest



Cancer Classifier

- MobileNetv2 (7 classes)
 - Default Batch normalization
 - 80% train / 20% validation
- Preprocessing
- Class Balancing
 - Oversampling of malignant class and class weights
- 15 Epochs
- Validation accuracy = 89%

Core ML 2-Stage Pipeline Evaluation

- Mole detector -> malignant classifier
- Critical oversight in malignancy detection (1k size sample)
 - Binary Accuracy: 84.10%
 - Cancer Recall: 11%
 - 89% of cancer cases were missed
 - Cancer Precision: 64%
 - 36% of cancer predictions were false alarms

Core ML 3-Stage Pipeline Evaluation

- Still critical, but better than before
- Confusion metrics (1k size sample):
 - Binary Accuracy: 46.70%
 - Cancer Recall: 30%
 - 70% of cancer cases were missed
 - Cancer Precision: 11%
 - 89% of cancer predictions were false alarms
 - Sensitivity over specificity for medical detection products

Future Improvements

- Improved accuracy
- Guide on taking/submitting ideal images
- Tell user next steps

Demo

Discussion

- Results
 - The task was more complex than we initially thought
 - Could possibly be viable for clinical use
 - Not prepared for real users
- Lessons Learned
 - Overfitting
 - Data augmentation
 - Data availability
 - Training techniques
- Limitations
 - Skin color
 - Malignancy type
 - On device computational limitations